




RETRACTORS



www.convertualinterlogic.com



+234 8033200998
+234 8180003461



enquiries@convertualinterlogic.com



1 Anu Crescent, Off Ado-Badore Road,
Ajah, Lagos, Nigeria



Fig. 2 30 x 10 mm
Fig. 1 30 x 13 mm

Fig. 1 26 x 10 mm
Fig. 2 34 x 13 mm

Farabeuf
J-19-001
12 cm / 4 3/4"
Set of 2



Fig. 2 27 x 16 mm
Fig. 1 26 x 16 mm

Fig. 1 23 x 16 mm
Fig. 2 30 x 16 mm

Farabeuf
J-19-002
15 cm / 6"
Set of 2



Fig. 2 21 x 14 mm
Fig. 1 17 x 14 mm

Fig. 1 28 x 16 mm
Fig. 2 32 x 16 mm

Collin
J-19-008
15 cm / 6"
Set of 2




11 x 4 mm

19 x 6 mm

Crile
J-19-010
11 cm / 4 1/4"



Mini-Roux
J-19-011
14 cm / 5 1/2"



Baby-Roux
J-19-012
13 cm / 5 1/4"



Roux

J-19-013 Set of 4
J-19-014 14 cm / 5 1/2"
J-19-015 15 cm / 6"
J-19-016 17 cm / 6 3/4"
J-19-017 18 cm / 7"

Fig. 2 22 x 15 mm
Fig. 1 19 x 15 mm

Fig. 1 19 x 15 mm
Fig. 2 22 x 15 mm



Parker
J-19-018
Set of 2




Fig. 2 29 x 25 mm
Fig. 1 26 x 25 mm

Fig. 1 24 x 25 mm
Fig. 2 27 x 25 mm

Parker

J-19-019

18 cm / 7"
Set of 2




Fig. 2 25 x 24 mm
Fig. 1 21 x 24 mm

Fig. 1 68 x 24 mm
Fig. 2 72 x 24 mm

Parker-Mott

J-19-020

16 cm / 6 1/4"
Set of 2



Fig. 2 25 x 24 mm
Fig. 1 21 x 24 mm

Fig. 1 40 x 15 mm
Fig. 2 45 x 15 mm

U.S. Army

J-19-024

21 cm / 8 1/4"
Set of 2




Fig. 2 47 x 13 mm
Fig. 1 41 x 13 mm

Fig. 1 37 x 25 mm
Fig. 2 42 x 25 mm

Mathieu

J-19-028

20 cm / 8"
Set of 2



36 mm

42 mm

Czerny

J-19-029

17.5 cm / 7"



Canny Ryoll

J-19-030 16 cm / 6 1/4"

J-19-031 20 cm / 8"

J-19-032 25 cm / 10"


J-19-033 30 cm / 12"



Goelet
J-19-034
19 cm / 7½"




40 x 40 mm



30 x 30 mm



Hibbs
J-19-037
22 cm / 8¾"
25 x 75 mm blade




J-19-035 115 x 27 mm
J-19-036 115 x 22 mm

Simon
28 cm / 11"



Taylor
J-19-038 18 cm / 7", 75 x 30 mm blade
J-19-039 20 cm / 8", 100 x 30 mm blade
Spinal retractor



Deaver
J-19-040
18 cm / 7"
19 mm

Deaver
J-19-041
21.5 cm / 8½"
22 mm

Deaver
J-19-042
23 cm / 9"
25 mm

Deaver
J-19-043
30 cm / 12"
25 mm

Deaver
J-19-044
33 cm / 13"
25 mm

Deaver
J-19-045
36 cm / 14"
25 mm




Deaver
J-19-046
30 cm / 12"
38 mm





Deaver
J-19-047
30 cm / 12"
50 mm



Deaver
J-19-048
30 cm / 12"
75 mm







Retractors

Micro Hooklets, Skin Hooks

Mikrohakchen, Hauthakchen


Micro Ganchitos, Ganchitos de piel





Desmarres

- J-19-086 8 mm
 - J-19-087 10 mm
 - J-19-088 12 mm
 - J-19-089 14 mm
 - J-19-090 16 mm
 - J-19-091 18 mm
- 16 cm / 6 1/4"



- J-19-092 12 mm
 - J-19-093 14 mm
 - J-19-094 16 mm
 - J-19-095 18 mm
- 16 cm / 6 1/4"



Langenbeck-Green
16 cm / 6 1/4"




Mini-Langenbeck


- J-19-098 10 x 6 mm
 - J-19-099 17 x 5 mm
 - J-19-100 20 x 6 mm
 - J-19-101 22 x 8 mm
- 16 cm / 6 1/4"




Cushing
20 cm / 8"



Cushing
24 cm / 9 1/2"







Ragnell-Davis
J-19-114
 14 cm / 5 1/2"



Senn-Green
 15 cm / 6"



Senn-Miller
 16 cm / 6 1/4"



Meyerding
 18 cm / 7"



J-19-119
 18 x 7 mm
 Fig. 1
 Blunt



J-19-120
 16 x 5 mm
 Fig. 2
 Toothed



J-19-121
 7 x 4 mm
 Fig. 3
 Toothed



J-19-122
 5 x 10 mm
 Fig. 4
 Sharp



J-19-123
 11 x 4 mm
 Fig. 5
 Blunt



J-19-124
 6 x 4 mm
 Fig. 6
 Saddle shaped



Gil-Vernet
 24 cm / 9 1/2"
 Flexible shaft



J-19-760
 8 mm



J-19-761
 11 mm



J-19-762
 15 mm





J-19-763
 17 mm



J-19-764
 20 mm



J-19-765
 23 mm





Murphy
19 cm / 7 1/2"
Sharp



J-19-615
2 prongs



J-19-616
3 prongs



J-19-617
4 prongs



J-19-618
6 prongs



Murphy
19 cm / 7 1/2"
Blunt



J-19-620
2 prongs



J-19-621
3 prongs



J-19-622
4 prongs



J-19-623
6 prongs



Kocher
22 cm / 8 1/2"



Sharp



J-19-145



J-19-146



J-19-147



J-19-148



J-19-149



J-19-150



Blunt



J-19-151



J-19-152



J-19-153



J-9-154



J-19-155



J-19-156




Koerte
J-19-157
40 x 30 mm
Sharp




Koerte
J-19-158
23.5 cm / 9 1/4"
40 x 30 mm
blunt




Israel
25 cm / 10"



J-19-162
70 x 70 mm




J-19-161
50 x 60 mm



J-19-160
45 x 50 mm



J-19-159
40 x 40 mm



Langenbeck
21 cm / 8 1/4"



J-19-163
30 x 11 mm



J-19-165
30 x 14 mm




J-19-167
30 x 16 mm



J-19-169
40 x 11 mm



J-19-171
50 x 11 mm



Langenbeck
22 cm / 8 1/2"



J-19-173
30 x 11 mm



J-19-175
30 x 14 mm




J-19-177
30 x 16 mm



J-19-179
40 x 11 mm





J-19-181
50 x 11 mm



Cushing
J-19-183
20 cm / 8"



Langenbeck
J-19-184 22 x 6 mm
J-19-185 44 x 13 mm



J-19-188	40 x 18 mm
J-19-189	60 x 20 mm
J-19-190	60 x 25 mm
J-19-191	75 x 30 mm
J-19-192	75 x 40 mm




Langenbeck
J-19-186
21 cm / 8¼"
85 x 15 mm



Langenbeck
J-19-187
21 cm / 8¼"
60 x 20 mm



Kocher
23 cm / 9"




J-19-193	25 x 6 mm
J-19-194	38 x 6 mm
J-19-195	35 x 11 mm
J-19-196	35 x 15 mm
J-19-197	40 x 11 mm
J-19-198	55 x 11 mm
J-19-199	70 x 14 mm
J-19-200	80 x 12 mm
J-19-201	80 x 16 mm

Kocher-Langenbeck
21.5 cm / 8 1/2"




Morris	
J-19-202	40 mm
J-19-203	50 mm
J-19-204	65 mm




Morris	
J-19-205	30 x 40 mm
J-19-206	45 x 50 mm
J-19-207	55 x 60 mm




Middeldorf



J-19-210
21.5 cm / 8 1/2"
14 x 17 mm



J-19-211
22 cm / 8 3/4"
20 x 22 mm



J-19-212
23.5 cm / 9 1/8"
26 x 30 mm



J-19-213
45 x 45 mm



J-19-214
45 x 50 mm




J-19-215
55 x 70 mm



Ollier
J-19-209
22.5 cm / 9"
36 x 30 mm



Ollier
J-19-208
22.5 cm / 9"
36 x 60 mm



Middeldorf
25.5 cm / 10"




J-19-216 28 x 20 mm
J-19-217 36 x 28 mm
J-19-218 44 x 38 mm



Richardson
24 cm / 9½"



Richardson
J-19-219
24 cm / 9½"
52 x 22 mm



Kelly
J-19-720 50 x 35 mm
J-19-720 65 x 50 mm
J-19-722 75 x 65 mm
26 cm / 10¼"



Fig. 2 49 x 38 mm
Fig. 1 28 x 20 mm


Fig. 1 36 x 28 mm
Fig. 2 71 x 43 mm

Richardson-Eastman
J-19-221 26 cm / 10½"
J-19-222 28 cm / 11"
J-19-800 Set of 2, Fig. 1 & Fig. 2







Tuffier
J-19-232
23 cm / 9"
64 x 60 mm



Koerte
J-19-233
28 cm / 11"
120 x 47 mm
Rigid



Heaney
J-19-234
27 cm / 10½"
82 x 28 mm



J-19-241
190 x 38 mm

J-19-242
190 x 57 mm



Mikulicz
J-19-237
26 cm / 10 1/4"
90 x 50 mm



Mikulicz
J-19-238 95 x 35 mm
J-19-239 120 x 50 mm
J-19-240 150 x 50 mm
26 cm / 10 1/4"



Kelly
27 cm / 10 1/2"



Meyerding
For deep surgery

J-19-245
25.5 cm / 10"
89 x 50 mm



J-19-244
23 cm / 9"
76 x 25 mm




J-19-243
21 cm / 8 1/4"
52 x 16 mm




Davidson
J-19-246 75 x 50 mm
J-19-626 75 x 90 mm
16 cm / 6 1/4"
Scapula retractor



Dartigues
J-19-247







Throkson
J-19-630 40 x 170 mm
J-19-631 50 x 170 mm
J-19-632 65 x 170 mm
33 cm / 13"




Throkson
J-19-640 40 x 170 mm
J-19-641 50 x 170 mm
J-19-642 65 x 170 mm
33 cm / 13"




Finsen
J-19-259 Sharp
J-19-260 Blunt
 7 cm / 2 3/4"




Finsen
J-19-262 Sharp
J-19-263 Blunt
 7 cm / 2 3/4"
 With adjusting screw



J-19-252
 8 cm / 3 1/4"



Buckley
J-19-251
 6 cm / 2 1/2"



J-19-253 7 cm / 2 3/4"
J-19-254 10 cm / 4"
J-19-255 Sharp 7 cm / 2 3/4"
J-19-256 Blunt 10 cm / 4"

Alm



J-19-458
 Sharp

J-19-458
 Blunt

10 cm / 4"




Jansen
 10 cm / 4"




J-19-264
 Sharp



J-19-265
 Blunt




Jansen-Wagner
J-19-257
 10 cm / 4"



Weitlaner
10.5 cm / 4 1/4"
2:3 Teeth




J-19-266
Sharp



J-19-267
Blunt




Weitlaner
3:4 Teeth



J-19-268 14 cm / 5 1/2"
J-19-270 16.5 cm / 6 1/2"
J-19-271 18 cm / 7"
J-19-272 20 cm / 8"
J-19-273 24.5 cm / 9 1/2"
Sharp




J-19-468
J-19-470
J-19-471
J-19-472
J-19-473
Blunt




West
J-19-274
13 cm / 5 1/4"
3:3 Teeth




Weitlaner-Wullstein
14 cm / 5 1/2"
3:3 Teeth




J-19-275
Sharp




J-19-276
Blunt




Weitlania
J-19-277
 14 cm / 5 1/2"



Mollison
J-19-278
 15 cm / 6"
 4-4 Teeth



Cloward
J-19-475
 13.5 cm / 5 1/4"
 Blunt




Milligan
J-19-279
 13 cm / 5"
 Blunt



Plester
J-19-280
 11 cm / 4 1/4"
 Blunt
 2 Teeth




Plester
J-19-281
 11 cm / 4 1/4"
 Blunt
 2 Teeth



Cloward
J-19-477
 13.5 cm / 5 1/4"
 Blunt



Perkins
J-19-282
13 cm / 5"



Perkins
J-19-283
13 cm / 5"



Plester
J-19-284
14 cm / 5 1/2"
Blunt



Norfolk-Norwich
J-19-285
22 cm / 8 1/2"
Blunt




Travers
J-19-286
20 cm / 8"



Adson-Baby / Cone
J-19-480
13 cm / 5 1/4"
Sharp




J-19-481
13 cm / 5 1/4"
Blunt



Adson-Baby / Cone
J-19-485
16 cm / 6 1/4"
Sharp



J-19-486
16 cm / 6 1/4"
Blunt



Adson / Cone
J-19-487
20 cm / 8"
Blunt




Adson

J-19-287	3 x 4	J-19-289
J-19-288	4 x 4	J-19-290
J-19-488	6 x 6	J-19-490
26.5 cm / 10 1/2"		26.5 cm / 10 1/2"
Sharp		Blunt



Adson Baby

J-19-291	14 cm / 5 1/2"	J-19-294
J-19-292	16 cm / 6 1/4"	J-19-295
J-19-293	20 cm / 8"	J-19-296
3-4 Teeth		3-4 Teeth
Sharp		Blunt



Harvey Jackson


J-19-297
28 cm / 11"
Sharp


Retractors

Wound Spreaders, Self-retaining

Wundspreizer

Separadores autoestaticos





Beckmann
J-19-298
 31 cm / 12 1/4"
 4-4 Teeth
 Sharp



Adson
J-19-299
 32 cm / 12 5/8"
 4.5 Teeth
 Blunt




Beckmann-Adson	
J-19-300	28 x 28 mm
J-19-301	44 x 28 mm
J-19-302	
J-19-303	
31 cm / 12 1/4"	31 cm / 12 1/4"
4-4 Teeth Sharp	4-4 Teeth Blunt

Retractors

Wound Spreaders, Self-retaining

Wundsperrizer

Separadores Autoestaticos




Beckmann-Eaton
J-19-304
32 cm / 12 5/8"
7:7 Teeth


Beckmann
J-19-305
31 cm / 12 1/4"



Gelpi
J-19-306 9 cm / 3½"
J-19-307 11.5 cm / 4½"
J-19-308 14 cm / 5½"
J-19-309 16 cm / 6¼"
J-19-310 18 cm / 7"
J-19-500 20 cm / 8"



Gelpi
J-19-311 9 cm / 3½"
J-19-312 11.5 cm / 4½"
J-19-313 14 cm / 5½"
J-19-314 16 cm / 6¼"
J-19-315 18 cm / 7"




Gelpi
J-19-316
 18 cm / 7"




Rigby
J-19-317
 18 cm / 7"



Mayo-Adams
J-19-318
 17 cm / 6¾"



Mayo-Adams
J-19-319
 17 cm / 6¾"
 With center blade




Meyerding

- J-19-320** 25 x 60 mm blade
- J-19-321** 25 x 65 mm blade
- J-19-322** 25 x 35 mm blade
- J-19-323** 22 x 60 mm blade




Darling

J-19-505
18 cm / 7"
(Poplitea)




Darling

J-19-510
18 cm / 7"
(Poplitea)




Spikel Shield

J-19-515 12.5 cm / 5"
J-19-516 18 cm / 7"
Left




Spikel Shield

J-19-520 12.5 cm / 5"
J-19-521 18 cm / 7"
Right



Collin, Baby
J-19-525
 20 cm / 8"
 Complete with 2 pairs of lateral blades
 (18 x 40 mm), (27 x 40 mm)




Oettingen
J-19-530
 Lateral Blades 38 x 80 mm
 Center blades 50 x 65 mm
 Blade opening 150 mm / 6"




Collin
J-19-328 With center blade
J-19-329 Without center blade
 25 cm / 10"
J-19-535 With center blade
J-19-536 Without center blade
 28 cm / 11"




Holzbach
J-19-337 Complete
J-19-338 Spreader only
J-19-339 59 x 80 mm lateral blades
J-19-340 74 x 82 mm center blade




Masson-Judd
J-19-326
Bladder retractor with
2 swivel blades 25 x 70 mm




Masson-Judd
J-19-327
Retractor with
2 swivel blades,
center blade 25 x 70 mm




Badgley
J-19-325
50 x 44 mm
16 cm / 6 1/4"
Spread



Frazier
J-19-324
8 cm / 3 1/4"
Spread
Laminectomy Retractor



Ricard
J-19-330 Complete set
J-19-331 Spreader only
J-19-332 55 x 70 mm center blade
J-19-333 50 x 58 mm lateral blades
28 cm / 11"
J-19-540 Complete set
J-19-541 Spreader only
J-19-542 80 x 90 mm center blade
J-19-543 60 x 80 mm lateral blades
30 cm / 12"




Rochard

J-19-341	48 x 90 mm
J-19-342	48 x 105 mm
J-19-343	48 x 120 mm
J-19-344	48 x 135 mm
J-19-345	48 x 155 mm

31 cm / 12 1/4"




Rochard
J-19-346
 19 cm / 7 1/2"
 Fixation device for
 J-19-341 - J-19-345



Doyen


J-19-347	48 x 90 mm
J-19-348	48 x 105 mm
J-19-349	48 x 120 mm
J-19-350	48 x 135 mm

28 cm / 11"




O'Sullivan-O'Connor


J-19-351	Adult
J-19-550	Baby




Gosset-Baby
J-19-353






Gosset
J-19-352




Gosset
J-19-354



Balfour-Baby
J-19-356 Single Bar
J-19-555 Double bar




Retractors
Abdominal Retractors
Bauchdeckenhalter
Separadores Abdominales




Balfour
J-19-359



Balfour
J-19-360 45 x 80 mm
J-19-361 70 x 85 mm




Balfour
J-19-362



Balfour	Opening
J-19-363	180 mm
J-19-560	200 mm
J-19-561	250 mm

With ratchet and detachable blades

Measurements of blades: Depth x Width




J-19-364
70 x 100 mm



J-19-365 75 x 45 mm
 J-19-366 65 x 45 mm
 J-19-367 82 x 70 mm



J-19-565
50 x 127 mm



Kirschner
 Standard, oval
J-19-370
 Frame only 300 x 240 mm

Kirschner
 Standard, oval
J-19-372
 Complete set, consisting 1 frame & 4 blades
 1 frame, 300 x 240 mm
 1 blade, 50 x 65 mm
 1 blade, 80 x 90 mm
 1 blade, 98 x 50 mm
 1 blade, 105 x 35 mm

Kirschner
 Standard, oval
J-19-750
 Complete set, consisting 1 frame & 7 blades
 1 frame, 300 x 240 mm
 1 blade, 40 x 55 mm
 1 blade, 50 x 65 mm
 1 blade, 60 x 85 mm
 1 blade, 78 x 65 mm
 1 blade, 80 x 90 mm
 1 blade, 98 x 50 mm
 1 blade, 105 x 35 mm

Measurements of blades: Depth x Width




J-19-374
 40 x 55 mm
 Standard




J-19-375
 50 x 65 mm
 Standard



J-19-376
 60 x 85 mm
 Standard



J-19-377
 78 x 65 mm
 For obese patients



J-19-378
 80 x 90 mm
 For obese patients




J-19-379
 98 x 50 mm
 For intestinal and bladder




J-19-380
 105 x 35 mm
 For intestinal and bladder

Abdominal Retractors for new-born and children
Extra blades are available only with set.



Kirschner
J-19-384
Complete set, consisting 1 frame & 4 blades
1 Frame, 240 x 240 mm
2 blades, 40 x 55 mm
2 blades, 50 x 65 mm


Kirschner
J-19-383
Frame only
240 x 240 mm



Kirschner
J-19-374
40 x 55 mm

Kirschner
J-19-375
50 x 65 mm

J-19-387



Denis-Browne
J-19-388
Complete set, consisting 1 frame & 4 blades
1 Frame, 150 x 175 mm
2 blades, 20 x 30 mm
2 blades, 20 x 40 mm

Denis-Browne
J-19-391
150 x 175 mm
Frame only




J-19-394
20 x 30 mm




J-19-395
20 x 40 mm

Extra blades are available only with set.



Franz
J-19-404
Complete set, consisting 1 frame & 8 blades
Folding frame, 300 x 300 mm
2 blades, 50 x 65 mm
2 blades, 60 x 85 mm
2 blades, 38 x 60 mm
2 blades, 45 x 80 mm



Franz
J-19-398
300 x 300 mm
Frame only
Folded



J-19-375
50 x 65 mm




J-19-401
38 x 60 mm



J-19-376
60 x 85 mm



J-19-402
45 x 80 mm



Scoville-Haverfield

J-19-576

Frame only, with joint
 Spreading 150 mm

J-19-575

Complete set of 13 pieces (hooks & blades)



Scoville

J-19-571

Frame only
 Spreading 150 mm

J-19-570

Complete set of 13 pieces (hooks & blades)

Blades on page 183

Blades for J-19-570 & J-19-575 on page 182



Scoville
J-19-580
 50 mm deep
 Blunt



Scoville
J-19-581
 70 mm deep
 Blunt



Scoville
J-19-585
 50 mm deep
 Sharp



Scoville
J-19-586
 70 mm deep
 Sharp



Scoville
J-19-590 60 x 25 mm
J-19-591 64 x 50 mm
J-19-592 88 x 25 mm
 Toothed



Taylor
J-19-600 48 x 50 mm
J-19-601 64 x 50 mm
J-19-602 64 x 64 mm
J-19-603 80 x 67 mm
 Sharp



Scoville
J-19-610 41 x 38 mm
J-19-611 67 x 44 mm
 Sharp



J-19-405
40 mm



J-19-408
40 mm



J-19-415
45 mm



J-19-409
45 mm



J-19-416
50 mm



J-19-410
50 mm



J-19-406
55 mm



J-19-411
55 mm



J-19-407
60 mm



J-19-412
60 mm



Cloward
16 cm / 6 1/4"
J-19-896
Complete Set



Cloward
25 cm / 10"
J-19-897
Complete Set